




# COOLPALL™ VERTOS ADVANCE



## OUR NEW BULK SHIPPING SYSTEM FEATURING PCM COOLANTS

CoolPall™ Vertos Advance is the NEW high performance bulk shipping system from Pelican BioThermal™. This bulk shipper features PCM coolants and has been designed to be flat packed, modular, and easily assembled. The CoolPall™ Vertos Advance maintains 2°C to 8°C for 120 hours of temperature controlled performance.

## WHY CHOOSE COOLPALL™ VERTOS ADVANCE?

### Ease & Speed of Use

- Superior performance foam and PCM technology reduces excursions and allows for hibernation to be achieved during transit, which pauses phase change within the shipper and can yield longer temperature compliant shipments
- Ergonomic design utilizes innovative interlocking components to ensure quick and easy assembly
- Modular components reduces inventory and allows for simple ordering of spare parts
- Reduced flat-packed height minimizes in-bound distribution costs and frees up valuable warehouse storage space
- Accepts full pallets of product for operational speed and simplicity

### Robust performance

- Provides strict temperature control of product within the demanding environment of global distribution
- Utilizes high density moulded materials for thermal and mechanical robustness while remaining lightweight
- Best in class payload to external volume ratio, lowering transport costs

## COOLPALL™ VERTOS ADVANCE

### DURATION

120+ Hours

### TEMPERATURE

2°C to 8°C



### VOLUME

ISO Pallet size,  
volume up to 1680L

# COOLPALL™ VERTOS ADVANCE PRODUCT SPECIFICATION

Product	Exterior Dimensions - mm	Interior Dimensions - mm	Payload Volume	Tare Weight - kg
ISOPAD	1355 x 1565 x 1595	1127 x 1296 x 1150	1680L	235

Peli BioThermal Profile, CoolPall™ Vertos Advance +2°C to +8°C



Interlocking high density construction



Inspection hatch

Minimizing the potential for excursions during customs inspection



FOR MORE INFORMATION ABOUT PELICAN BIOTHERMAL PRODUCTS, PLEASE CONTACT:

INFO@PELICANBIOTHERMAL.COM T: +1 763 412 4800 PELICANBIOTHERMAL.COM

ISO17025, ISO9001-2008 Accredited

Golden Hour™, TIC™, Crêdo™, Crêdo Cube™, Crêdo ProEnvision™, Crêdo Xtreme™, Crêdo Duracube™, Crêdo ProMed™, Sherpa Systems™, CoolGuard™ and CoolPall™ are registered trademarks of Pelican BioThermal LLC., its subsidiaries and/or affiliates.



©2020 PELICAN BioThermal all rights reserved.